

Evaluation by real-time PCR:

FastGene® PCR adhesive clear seal (FG-93AC) and FastGene® qPCR pressure clear seal (FG-95PC).

Purpose: Test whether FastGene® PCR adhesive clear seal is suitable for real-time PCR; Test was carried out in comparison with the commercially available competitor Fs qPCR pressure clear seal.

Evaluated products: FastGene® PCR adhesive clear seal (FG-93AC) and FastGene® qPCR pressure clear seal (FG-95PC)

Cat. No. FG-93AC; FG-95PC

PCR conditions

qPCR instrument: StepOnePlus Life Technologies (ABI)
Plate: Competitor Fs 96-well plate
Evaluated Seals: 1) FastGene® PCR adhesive clear seal (FG-93AC)
 2) FastGene® qPCR pressure clear seal (FG-95PC)
Comparison product: 3) Competitor Fs qPCR pressure clear seals
qPCR reagent: KAPA SYBR Fast qPCR Kit
Template DNA: Roche Human Genomic DNA (# 1169112001) 5.0 ng / µl
Primer: Act-F1, Act-R1 (10 µM) < β-actin: 294 bp amplicon

DNA conc. [µg/µL]	1	2	3	4	5	6	7	8	9	10	11	12
A	●	●	●	●	●	●	●	●	●	●	●	●
B	●	●	●	●	●	●	●	●	●	●	●	●
C	●	●	●	●	●	●	●	●	●	●	●	●
D	●	●	●	●	●	●	●	●	●	●	●	●
E	●	●	●	●	●	●	●	●	●	●	●	●
F	●	●	●	●	●	●	●	●	●	●	●	●
G	●	●	●	●	●	●	●	●	●	●	●	●
H	●	●	●	●	●	●	●	●	●	●	●	●

Act-F1: TCACCCACACTGTGCCCATCTACGA
Act-R1: CAGCGGAACCGCTCATTGCCAATGG

Reaction components

2xMaster Mix 10 µl
 H₂O 7.6 µl
 Primer F (10µM) 0.2 µl
 Primer R (10µM) 0.2 µl
 Template DNA 2.0 µl
 total 20 µl

Cycle program (StepOne default fast cycle)

Initial Denaturation: 95 °C 20 sec
 Denaturation: 95 °C 3 sec
 Anneal/Extension: 60 °C 30 sec (detection) 40 cycle
 Melting Curve: 95 °C 15 sec
 60 °C 1 min
 95 °C 15 sec (+0.3 °C step detection)
 Cycle simulation: 20 µl

Evaluation of appearance



Evaluated products

Comparison product

1) FastGene® PCR adhesive clear seal (FG-93AC)

2) FastGene® qPCR pressure clear seal (FG-95PC)


3) Competitor Fs qPCR pressure clear seal



Compared to the both pressure seals, irregularities at the adhesive clear were observed on the surface of the seal. This irregularity is presumably caused by the adhesive material (glue).


Evaluation of amplification curves

1) FastGene® PCR adhesive clear seal (FG-93AC)




Conc. [ug/uL]	Ct	Ct SD
5000	23.3	0.09
1000	25.82	0.08
200	28.31	0.07
40	30.96	0.11

2) FastGene® pressure clear seal (FG-95PC)



Conc. [ug/uL]	Ct	Ct SD
5000	23.38	0.07
1000	25.88	0.07
200	28.46	0.09
40	31.04	0.18

3) Competitor Fs pressure clear seal




Conc. [ug/uL]	Ct	Ct SD
5000	23.3	0.09
1000	25.81	0.12
200	28.31	0.11
40	30.89	0.09


The qPCR amplification curves are in the same range for all three seals. No significant difference could be observed between the seals. The comparison of the Ct and Ct SD show for all three seals very similar values.

Evaluation of the melting curve


1) FastGene® PCR adhesive clear seal (FG-93AC)



2) FastGene® pressure clear seal (FG-95PC)




3) Competitor Fs pressure clear seal




The results for the melting curve are almost the same.

Evaluation of the calibration curve


1) FastGene® PCR adhesive clear seal (FG-93AC)



2) FastGene® pressure clear seal (FG-95PC)



3) Competitor Fs pressure clear seal



The calibration curve of the three different sheets show no difference in regards to slope, efficiency or variation

Conclusion

According to the results, FastGene® PCR adhesive clear seal (FG-93AC) showed similar performance compared to the other two qPCR pressure seals. The seal can therefore be used without restrictions in our qPCR experiments.