




Application

Quantitative analysis of gene expression in human bone marrow-derived mesenchymal stem cells (primer check)

Product Name	FastGene™ RNA Premium Kit (FG-81050, FG-81250)
Manufacturer	NIPPON Genetics EUROPE
Product Name	Transcriptor Universal cDNA Master (05893151001)
Manufacturer	Roche Diagnostics
Product Name	KAPA SYBR® FAST qPCR Kits (KK4610, KK4611)
Manufacturer	Kapa Biosystems, Inc.

The data below was published through the courtesy of Mr. Mitsuru Furuya, Advanced Polymeric Nanostructured Materials Engineering, Department of Advanced Science and Technology, Toyota Technological Institute, Toyota School Foundation, Japan.

Total support with our products was possible!



Purposes of these experiments

1) In preparation for planned experiments in which hMSC (human mesenchymal stem cells) would be induced to differentiate into cartilage and bone, we checked primers targeting the following 10 genes using hMSC-BM (human bone marrow-derived mesenchymal stem cells).

Genes: ACTB, GAPDH, RUNX2 (2 types), ALPL, Collagen Type I, SP7/Osterix, BGLAP (2 types), SOX9, Aggrecan, and Collagen Type II (11 types)

2) We verified whether 2.5 µg of RNA, which is a recommended maximum amount of input RNA for reverse transcription reaction, can be used for qPCR analysis (Sample A).

Experimental conditions

■ RNA extraction : FastGene™ RNA Premium Kit

- Starting sample: hMSC-BM (human bone marrow-derived mesenchymal stem cells) 1.8×10^6 cells
- The amount of lysis buffer RL: 600 µL
- Pretreatment of the sample: Lysis buffer RL (600 µL) was added onto the cells, and the cells were not homogenized but lysed by vortexing alone.
- Elution buffer RE: 50 µL
- Measurement of yield and purity: NanoDrop One (Thermo Fisher Scientific, Inc.)

■ Reverse transcription : Transcriptor Universal cDNA Master

- The amount of input RNA: Total RNA 1.0 µg (Sample B) or 2.5 µg (Sample A)
- The type of primers: Random hexamer primers
- The composition of the reaction mixture:

RNA	1.0 µg or 2.5 µg
Transcriptor Universal Reaction Buffer	4 µL
Transcriptor Universal Reverse Transcriptase	1 µL
Distilled Water	up to 20 µL
total	20 µL

- Reaction program:
- | | |
|-------------------------------|---------------|
| Annealing | 25°C · 5 min |
| ↓ | |
| Enzyme reaction | 55°C · 10 min |
| ↓ | |
| Denaturation and inactivation | 85°C · 5 min |



■ Quantitative analysis of gene expression : KAPA SYBR® FAST qPCR Kit

- The amount of input cDNA: 2 μ L (per 10- μ L reaction) of the reverse transcription reaction solution diluted by 10-fold with DEPC-treated water
- Real-time machine: LightCycler® 96
- The composition of the reaction mixture:

cDNA	2 μ L
KAPA SYBR® FAST qPCR Master Mix (2X)	5 μ L
10 μ M forward primer	0.2 μ L
10 μ M reverse primer	0.2 μ L
PCR-grade water	Up to 10 μ L
total	10 μL

• Reaction program:

Initial denaturation	95°C · 3 min
↓	
Denaturation	95°C · 10 sec
Annealing	60°C · 20 sec
Extension reaction	72°C · 1 sec
↓	
Melting	
↓	
Cooling	40°C · 10 sec

Results

1. RNA extraction


The volume of elution buffer (μ L)	Concentration (ng/ μ L)	A260/A280	A260/A230	A260	A280
50	474.82	2.07	2.23	11.87	5.74

2. Quantitative analysis of gene expression

Sample A (RNA 2.5 μ g)

As shown below, all the primers detected respective targets.

● Amplification curve



● Cq values and positive/negative judgements


Primers	No.	Cq value	Result	Primers	No.	Cq value	Result	Primers	No.	Cq value	Result	Primers	No.	Cq value	Result
GAPDH	1	17.85	Positive	SP7	1	33.48	Positive	Col2 9/4_2	1	27.30	Positive	Col2 A	1	28.41	Positive
	2	17.89	Positive		2	32.48	Positive		2	27.30	Positive		2	28.39	Positive
ACTB	1	15.28	Positive	Bglap1	1	28.53	Positive	Col2 9/4_3	1	31.76	Positive	Col2 B	1	32.38	Positive
	2	15.27	Positive		2	28.66	Positive		2	31.53	Positive		2	34.27	Positive
RUNX2_1	1	22.98	Positive	Bglap3	1	30.24	Positive	Col2 9/4_4	1	31.32	Positive	Col2 C	1	33.62	Positive
	2	22.98	Positive		2	30.73	Positive		2	30.95	Positive		2	33.17	Positive
RUNX2_2	1	23.41	Positive	SOX9	1	22.53	Positive	Col2 8/30_1	1	39.03	Positive	Col2 D	1	33.48	Positive
	2	23.39	Positive		2	22.54	Positive		2	31.74	Positive		2	33.41	Positive
ALPL	1	25.59	Positive	Aggrecan	1	37.41	Positive	Col2 8/30_2	1	31.43	Positive				
	2	25.63	Positive		2	35.06	Positive		2	32.34	Positive				
Col1	1	14.26	Positive	Col2 9/4_1	1	31.00	Positive	Col2 8/30_3	1	35.11	Positive				
	2	14.30	Positive		2	30.48	Positive		2	37.93	Positive				



Sample B (RNA 1.0 µg)

As shown below, most of the primers detected respective targets.

● Amplification curve



● Cq values and positive/negative judgements

Primers	No.	Cq value	Result	Primers	No.	Cq value	Result	Primers	No.	Cq value	Result	Primers	No.	Cq value	Result
GAPDH	1	18.66	Positive	SP7	1	34.60	Positive	Col2 9/4_2	1	28.12	Positive	Col2 A	1	29.66	Positive
	2	18.67	Positive		2	-	Negative		2	28.17	Positive		2	30.15	Positive
ACTB	1	16.24	Positive	Bglap1	1	30.51	Positive	Col2 9/4_3	1	32.22	Positive	Col2 B	1	34.47	Positive
	2	16.22	Positive		2	29.61	Positive		2	32.47	Positive		2	34.67	Positive
RUNX2_1	1	24.15	Positive	Bglap3	1	30.95	Positive	Col2 9/4_4	1	33.04	Positive	Col2 C	1	36.38	Positive
	2	24.21	Positive		2	31.28	Positive		2	31.69	Positive		2	34.41	Positive
RUNX2_2	1	24.60	Positive	SOX9	1	23.90	Positive	Col2 8/30_1	1	32.56	Positive	Col2 D	1	34.11	Positive
	2	24.61	Positive		2	23.81	Positive		2	34.34	Positive		2	33.22	Positive
ALPL	1	26.45	Positive	Aggrecan	1	34.99	Positive	Col2 8/30_2	1	31.92	Positive				
	2	26.61	Positive		2	-	Negative		2	31.29	Positive				
Col1	1	15.39	Positive	Col2 9/4_1	1	30.70	Positive	Col2 8/30_3	1	34.37	Positive				
	2	15.37	Positive		2	30.85	Positive		2	33.96	Positive				

General evaluation

- LightCycler® 96 system (high performance, simple operation)
 - · Because its melt curve analysis does not take time, the experimental time can be reduced, which is good.
 - Because of its simple operation using a touchscreen, the machine is easy for beginners to use.
 - Its software is easy to use, for example it is compatible with Excel.
- FastGene™ RNA Premium Kit (low price, high quality)
 - · To be honest, I felt it burdensome to eliminate gDNA after RNA was eluted, instead of doing it on column. However, because you do it carefully, even when you are not accustomed to it, there is a sense of security that gDNA can be sufficiently eliminated.
 - Although there are multiple columns, they were color-coded in the protocol and easy to identify.
 - Reaction time is much reduced than the previous products, and it became easy to process many samples.
 - I am satisfied because its quality is equal to or better than Qiagen's RNeasy Mini Kit and its price is cheaper.
- Transcriptor Universal cDNA Master (convenience, high performance)
 - · cDNA can be easily synthesized by mixing 2 premixed reagents and RNA, and is high-quality. I am very satisfied. Because all of our laboratory members are beginners of PCR, I am thankful that it can be done by simple operation.
- KAPA SYBR® FAST qPCR Kit (low price, high sensitivity)
 - · Clearer calibration curve can be drawn than Cells to Ct kit



Customers' comments

When we purchased LightCycler®, we consulted Nippon Genetics, which is selling reagents from Roche, and they recommended these products because of their low price and high quality. Although the RNA purification step is burdensome, I am satisfied because RNA has been sufficiently purified in our laboratory which has many beginners.