

# Comparing Midori Green Advance to Competitor RS

Cat. No. MG04

## Introduction

Next generation nucleic acid dyes are a safe alternative for the carcinogenic Ethidium Bromide. Nonetheless, the dyes have different characteristics. Midori Green Advance is a dye from the latest generation being a further development of the original Midori Green, with better signal to background ratio than many competitors.

## Methods

Two gels were prepared and loaded with equal amounts of DNA. 35 ml of 1% Agarose were prepared with 1 x TAE Buffer (IDsol, Cat. No. #ID1521). Competitor RS was diluted 50000 x, while Midori Green Advance was diluted 25000 x.

Electrophoresis was performed using the Mupid ONE electrophoresis system (Advance Inc., Cat. No. #MU2). Two images were taken after 20 minutes and 40 minutes at 100V, respectively.

## Samples/lane:

- I. FastGene (Cat. No. MWP100)
- II.PCR fragment of 80 bp
- III.PCR fragment of 80 bp
- IV.Kapa Express ladder

## Conclusions

1. **Stronger signals**
2. **Better signal to background ratio**
3. **No destaining of small fragments**

## Results

The comparison of the gels showed that **Competitor RS has a substantially higher background** when compared to Midori Green Advance (Fig.1, white arrows). The DNA signal was therefore harder to detect.

After 40 minutes electrophoresis, **the staining of the DNA was more stable with Midori Green Advance than with Competitor RS** as seen after 40 minutes electrophoresis (Fig.1, red and green arrows). **The small DNA fragments disappeared when stained with Competitor RS (Fig.1, zoomed image)**. Therefore, the staining of small DNA fragments using Competitor RS has to be seen critically, as false negative results could lead to wrong scientific conclusions.

**Midori Green Advance is the superior DNA dye**




Figure 1: : Images of the gels after 20 minutes and 40 minutes electrophoresis at 100V.