



Application

RNA isolation of Neuro 2a cells (mouse neuroblastoma) Analysis of purified RNA by Real Time PCR (qPCR)

Product name

FastGene® RNA Basic Kit (FG-80050, FG-80250)

Manufacturer

Nippon Genetics Europe

The following data was published due to the kindness of the Department of Peripheral Nervous System Research, National Institute of Neuroscience, National Center of Neurology and Psychiatry, Japan

Method

Mouse derived neuroblastoma Neuro2a cells were divided in 6 parts (per tube 1 Mio cells). Three tubes were used for purification by the Qiagen RNeasy kit and the other three tubes were used with the FastGene® RNA Basic Kit.

1 µg of RNA was reverse transcribed in a 20 µl reaction system using ReverTra Ace enzyme of TOYOBO. The reverse transcription reaction sample was diluted 20-fold with MQ (Millipore) water and 1 µl was subjected to qPCR. For the PCR reaction the THUNDERBIRD SYBR qPCR Mix from TOYOBO Co., Ltd. was used. Actin and Gene A were used as target genes. All experiments were performed using the 7300 Real-Time PCR System from Applied Biosystems.

Condition

1. RNA purification

Sample type: Neuro2a cells (mouse derived neuroblastoma)
 Sample amount: 1.0×10⁶ cells / prep
 RNA purification kit: ①Company Q Kit ②FastGene® RNA Basic Kit (each n=3)
 DNase I treatment: Neither ① nor ②
 RNA elution volume: 30 µl for both kits

2. RT reaction

Input RNA amount: 1 µg (per 20 µl per reaction)
 Reverse transcriptase: ReverTra Ace (TOYOBO)

Reaction set up:

total RNA	1µg
Random Primers (25 pmoles/µl)	1µl
5×Buffer	4µl
10 mM dNTPs	2µl
ReverTra AceR (100units/µl)	1µl
(Total Volume)	20µl)

Set up:

Annealing 30°C•10min
 ↓
 Enzyme reaction 42°C•30min
 ↓
 Denaturation 99°C•5min


3. Real time PCR (qPCR)

cDNA volume: 1 µl of reverse transcription reaction solution diluted 20-fold (20 µl per 1 reaction)
 qPCR Reagent: THUNDERBIRD SYBR qPCR Mix (TOYOBO)
 qPCR Instrument: Applied Biosystems 7300 Real-Time PCR System

• n=3
 • DNase I not performed
 • 30µl elution

• RNA sample 1 µg/ 20 µl reaction

• 20-fold dilution cDNA template 1 µl/ 20 µl reaction
 • n=3



Result

Yield and purity

		yield		purity	
RNA sample		(µg)	Mean (µg)	A260/280	A260/230
Company Q	1	17.6625	18.7547	2.05	2.31
	2	18.4535		2.05	2.47
	3	20.1480		2.04	2.50
FastGene®	1	18.9745	19.8202	2.08	2.33
	2	21.4340		2.06	2.49
	3	19.0520		2.07	2.42

NanoDrop
(Thermo Scientific™ Ltd)

There is just a little difference between both kits in regard to yield and purity.

qPCR result: **actin**

Sample	Ct	Mean	Mean Ct	StdDev Ct	P value
company Q	1	16.1606	16.1369	16.3551	0.29728
		16.2475			
		16.0027			
	2	16.2128	16.2347		
		16.3102			
		16.1810			
	3	16.7193	16.6937		
		16.5803			
		16.7815			
FastGene®	1	15.9163	15.8943	15.8970	0.05795
		15.8641			
		15.9024			
	2	15.9028	15.9563		
		15.8711			
		16.0950			
	3	15.9139	15.8405		
		15.8597			
		15.7479			

qPCR result 2: **Gene A**

Sample	Ct	Mean	Mean Ct	StdDev Ct	P value
company Q	1	21.7528	21.7298	21.9913	0.34464
		21.5308			
		21.9058			
	2	22.3860	22.3819		
		22.5230			
		22.2366			
	3	21.5057	21.8623		
		22.2444			
		21.8369			
FastGene®	1	21.4406	21.4551	21.5940	0.13235
		21.5040			
		21.4208			
	2	21.2879	21.6081		
		21.8181			
		21.7182			
	3	21.9397	21.7187		
		21.5256			
		21.6908			

In comparison with the company Q kit, the result of qPCR show little difference ($P > 0.05$).



Customers comment

We have used Q products so far, but the FastGene® RNA Basic Kit from Nippon Genetics is very convenient to use and yield and purity are even a little bit better.