



Application

## Comparative evaluation of FastGene® Scriptase Basic and other company's kit by RT-PCR

Category

Reverse Transcriptase

Product

FastGene® Scriptase Basic (LS52)

Manufacturer

NIPPON Genetics EUROPE

The following data was published due to the kindness of Mr. Tamura, a member of the Japanese public research institute.

### Overview



#### FastGene® Scriptase Basic (LS52)

In order to evaluate the performance of FastGene® Scriptase Basic, reverse transcription was performed on the RNA extracted from HeLa cells with this product and other company's reverse transcriptase.

The obtained cDNA was amplified by PCR for the target gene and confirmed by electrophoresis. The performance of FastGene® Scriptase Basic was evaluated by comparing the band signals.

### Result

- RNA extraction

Sample: HeLa cells were cultured on a 35 mm dish and used at 70-80% confluency.

RNA extraction: Ribozol (MS Technosystems) (N580-100ML)

Absorbance measurement: U-3900 (HITACHI) (U-3900)

Result: A260/280 = 1.56, yield = 9.5 µg

- Reverse transcription reaction

Total RNA was diluted to 1, 10, 100 ng/µL. 1 µL of each dilution was added as Template RNA. ※ 0 ng is 1 µL H<sub>2</sub>O

- Primer

Gene specific primer (GAPDH): CTCTTCCTCTTGTGCTCTTGC

#### FastGene® Scriptase Basic

Template RNA (0, 1, 10, 100 ng/µL)	: 1 µL
Gene specific primer (GAPDH)	: 2 µL (2 pmol)

↓  
42 °C 5 min

↓  
add

10X FastGene® Scriptase buffer	: 2 µL
dNTP Mixture (2 mM each)	: 2 µL
FastGene® Scriptase basic	: 1 µL
Sterile water (RNase free)	: 12 µL

---

total : 20 µL

↓  
42 °C 60 min  
↓  
90 °C 5 min  
↓  
keep at 22 °C

#### company T, RT rection kit

Template RNA (0, 1, 10, 100 ng/µL)	: 1 µL
Gene specific primer (GAPDH)	: 2 µL (2 pmol)
2 mM dNTP	: 5 µL
RNase free dH <sub>2</sub> O	: 2 µL

---

total : 10 µL

↓  
65 °C 5 min

↓  
4 °C 5 min

↓  
add

5X buffer	: 4 µL
Reverse transcriptase	: 1 µL
RNase free dH <sub>2</sub> O	: 5 µL

↓  
42 °C 50 min  
↓  
70 °C for 15 min  
↓  
keep at 22 °C

Reference: Template total RNA recommended amount: 1 ng-5 µg

Reference: Template total RNA recommended amount: 5 µg or less



### PCR conditions

#### ● Reaction composition

10 μM Forward Primer	: 2.5 μL
10 μM Reverse Primer	: 2.5 μL
2 mM dNTPs (TOYOBO)	: 5 μL
5X Phusion HF Buffer	: 10 μL
cDNA	: 2 μL
sterile water	: 27.5 μL
Phusion DNA Polymerase (NEB)	: 0.5 μL
total	: 50 μL

#### ● Program

95 °C	2 min	
↓		
95 °C	10 sec	} 30 cycles
65 °C	30 sec	
72 °C	30 sec	
↓		
72 °C	3 min	
↓		
hold	22 °C	


#### ● Primers for-GAPDH

Forward Primer : CCACAGTCCATGCCATCAC  
 Reverse Primer : CCATGAGGTCCACCACCC  
 Amplification product size : 500bp

### Electrophoresis condition

Agarose gel	: 1% agarose gel
Sample	: PCR reaction solution (15 μL)
Running buffer	: 1xTAE
Migration conditions	: 100 V, 30 min
Marker	: 7 μL (Gene Ladder Wide 1) (Nippon Gene) (Code No. 313-06961)

### Result



We were able to detect PCR bands by using Scriptase Basic and a small amount of total RNA.



#### Customer's comment

Scriptase Basic is relatively inexpensive, so it is cost effective and satisfactory.  
 The tube labels are color-coded, which makes them easy to use.