



Technical Data

FastGene® 50bp DNA Ladder migration test

Evaluation product

FastGene® 50bp DNA Ladder (MWD50)

Evaluation method


Evaluation of electrophoretic photographs under different conditions of recommended use of FastGene® 50bp DNA Ladder.

* It is important to detect the low molecular region (50 - 250 bp)

Purpose

The recommended migration condition of FastGene® 50bp DNA Ladder (MWD50) was evaluated

Appliances and reagents



Fas-Digi Dark box body only
(Cat No. FAS-DGMU)

Fas-Digi exclusive digital camera
(Pentax MX1)

Blue/Green LED Illuminator (500nm)
(Cat No. FG-08)

Blue LED Illuminator (470nm)
(Cat No. FG-06)

UV transilluminator medium wavelength (302nm)
(Cat No. FG-300)




FastGene®
Agarose
(Cat No. AG01/AG02)



Staining reagent
Midori Green Direct
(Cat No. MG06)

Safe Blue Electrophoresis
Full System



FastGene®
50bp DNA Ladder
(Cat No. MWD50)

Result

FastGene® 50bp DNA (MWD50)

Electrophoretic picture	TAE		TBE	
	ng	bp	ng	bp
	40	1500	40	1500
	70	1200	70	1200
	33	1000	33	1000
	27	800	27	800
	23	700	23	700
	20	600	20	600
	67	500	67	500
	23	450	23	450
	27	400	27	400
	18	350	18	350
	20	300	20	300
	25	250	25	250
	67	200	67	200
	30	150	30	150
	27	100	27	100
30	50	30	50	
Buffer	TAE		TBE	
Agarose concentration	3.0%		3.0%	
Agarose capacity	12.5mL		12.5mL	
Gel size	H5.9cm×W5.1cm		H5.9cm×W5.1cm	
Electrophoresis time	60min		60min	
Voltage	100V		100V	
Staining reagent	Midori Green Direct		Midori Green Direct	
Agarose	FastGene® Agarose (Cat No. AG01/AG02)		FastGene® Agarose (Cat No. AG01/AG02)	

Recommended migration condition

TAE : Agarose concentration 3.0% Electrophoresis time 60min Staining reagent MGD


TBE : Agarose concentration 3.0% Electrophoresis time 60min Staining reagent MGD

Summary

TAE or TBE whichever is used, it is recommended to use an agarose concentration of 3% and an electrophoresis time of about 60 minutes.



Reference Data



* If the migration time is short, the low molecular area looks smeared.

Information

Tris / acetic acid / EDTA buffer (TAE)

<Characteristic>

- Most common
- Suitable for linear, supercoiled DNA isolation
- Approximately 10% faster migration than TBE
- pH variation tends to occur during electrophoresis since buffering ability is low. There is a need in buffer exchange or buffer circulation under high current conditions.
- x 50 available

Tris / boric acid / EDTA buffer (TBE)

<Characteristic>

- It has a strong buffering power, can suppress pH variation and heat generation, suitable for electrophoresis under high voltage conditions
- The migration time becomes longer than TAE
- May have consequences in downstream applications

