

Technical Data Sheet

FastGene® RNA Basic / Premium Kit DNase I Protocol Treatment

- **Goal:** Evaluation of the DNase I treatment after elution recommended in the FastGene® RNA Premium, and comparison to an on-column DNase I treatment.
- **Materials:**

| | | |
|---------------------------|-----------|-----------------|
| FastGene® RNA Basic kit | 6 Preps | Cat.No.FG-80006 |
| | 50 Preps | Cat.No.FG-80050 |
| | 250 Preps | Cat.No.FG-80250 |
| FastGene® RNA Premium kit | 6 Preps | Cat.No.FG-81006 |
| | 50 Preps | Cat.No.FG-81050 |
| | 250 Preps | Cat.No.FG-81250 |

Background

A DNase I treatment is not obligatory when using RNA purified with a silica membrane. However, for very DNA-sensitive downstream application, one of the following methods can be performed:

1. DNase I Treatment after elution: this is the standard protocol of the FastGene® RNA Premium kit.
2. DNase I Treatment on column: optional protocol available


<DNase I treatment after elution>

In this method, impurities, such as salts, are removed during the RNA purification. This enables an enzymatic reaction in ideal conditions with high DNA removal efficiency.

<DNase I treatment on-column>

This method is widely used, due to its convenience. Here, the DNase I treatment is performed on column after binding the RNA to the column. This is, however, in a high-salt concentration environment, which affects the DNase I treatment efficiency. In order to avoid this, the column must be washed with adequate washing buffer. Failed removal of salts will result in a lower enzymatic activity and too low concentration of salts will cause the release of the RNA, resulting in reduced yields.

The following DNA removal efficiency is expected as follows:



Here, we present results of four different approaches to remove DNA and describe the observed efficiencies.

Experimental Conditions

Sample: Jurkat cell line
 5× 10⁵ cells /prep
 n=3

- Evaluated points:
- 1 Yield
 - 2 RIN Score
 - 3 Residual genomic DNA rate


Test conditions

| Kit name | DNase Treatment |
|---------------------------|--------------------------------|
| FastGene® RNA Basic kit | ① None |
| | ② On -column |
| FastGene® RNA Premium kit | ③ Filtration and on-column |
| | ④ Filtration and after elution |
| Competitor Q | ① None |
| | ② On column |
| | ③ Filtration and on-column |

Results

Yield

| Product | DNase I Treatment | Results | |
|---------------------------|--------------------------------|---------------|-------------|
| | | Average Yield | Stand. Dev. |
| FastGene® RNA Basic kit | ① None | 3.81 | 0.07 |
| | ② On -column | 3.15 | 0.22 |
| FastGene® RNA Premium kit | ③ Filtration and on-column | 2.56 | 0.47 |
| | ④ Filtration and after elution | 2.27 | 0.32 |
| Competitor Q | ① None | 4.26 | 0.18 |
| | ② On column | 2.97 | 0.18 |
| | ③ Filtration and on-column | 2.85 | 0.24 |



All kits showed similar results under similar conditions.


Yields were in the following order:

No DNase I treatment > On-column treatment > Filtration of sample before on-column treatment > Filtration before DNase I treatment after elution

Reason for this difference in the yield measurement could be the presence of residual genomic DNA.

RIN Score

| Product | DNase I Treatment | Results | |
|---------------------------|--------------------------------|---------------|-------------|
| | | Average Yield | Stand. Dev. |
| FastGene® RNA Basic kit | ① None | 10 | 0 |
| | ② On -column | 10 | 0 |
| FastGene® RNA Premium kit | ③ Filtration and on-column | 10 | 0.2 |
| | ④ Filtration and after elution | 10 | 0 |
| Competitor Q | ① None | 9.8 | 0.2 |
| | ② On column | 9.8 | 0.3 |
| | ③ Filtration and on-column | 9.5 | 0 |




RIN score showed no difference under any condition.

Residual genomic DNA rate

Calculation of residual genomic DNA using absorbance and real-time PCR results:

$$\text{Residual genomic DNA\%} = \frac{\text{Amount of residual genomic DNA measured by qPCR [ng]}}{\text{Amount of RNA measured by absorbance [ng]}} \times 100$$

| Product | DNase I Treatment | Results | |
|---------------------------|--------------------------------|-------------------------|-------------------------|
| | | Average Yield | Stand. Dev. |
| FastGene® RNA Basic kit | ① None | 1.16 × 10 | 1.65 |
| | ② On -column | 1.56 | 8.07 × 10 ⁻² |
| FastGene® RNA Premium kit | ① Filtration and on-column | 1.18 × 10 ⁻² | 3.00 × 10 ⁻³ |
| | ④ Filtration and after elution | 9.74 × 10 ⁻⁴ | 5.41 × 10 ⁻⁴ |
| Competitor Q | ① None | 1.82 × 10 | 2.98 |
| | ② On column | 1.63 × 10 ⁻² | 1.47 × 10 ⁻³ |
| | ③ Filtration and on-column | 1.33 × 10 ⁻² | 9.13 × 10 ⁻⁴ |



We can confirm the assumption stated by our company in the beginning: A DNase I treatment after elution showed the lowest amount of residual genomic DNA with a higher reproducibility, when compared to the other tested conditions

• Summary

Based on this experiment, the following results were obtained

• Yield:

① Without DNase I treatment > ② On-column DNase I treatment > ③ pre-filtration + on-column DNase I treatment > ④ pre-filtration + DNase I treatment after elution

• RIN Score:

The same tendency under any condition

• Genomic DNA removal efficiency:

① Without DNase I treatment < ② On-column DNase I treatment < ③ pre-filtration + on-column DNase I treatment < ④ pre-filtration + DNase I treatment after elution

Based on these results, the "DNase I treatment after elution" is recommended for the FastGene® RNA Premium Kit, for being the most effective for genomic DNA removal.



NIPPON Genetics EUROPE
Binsfelder Straße 77, D-52351 Dueren Germany

<http://www.nippongenetics.eu>

TEL +49 2421554960 FAX +49 24215549611 info@nippongenetics.eu