

Quick Start High Resolution Quantitative Cartridge Kit

A. Specifications

Specifications	Description
Sample Volume Requirement	≥ 2 µl
DNA Sizing Range	20-1500 bp
DNA Sizing Precision*	2% CV
DNA Fragment Concentration Range	**0.2 ng/ µl to 50 ng/ µl
DNA Quantification Precision	10% CV

* Determined by the DNA ladder (15-622 DNA Size Marker, C109200) as sample

** Best Quantitative Range: 0.5-2 ng/µl

B. Kit Components and storage conditions

Contents	Storage Condition
High Resolution Cartridges (S1)	≥4°C
Quantitative Marker (20bp&1500bp, C109109-500AQ)	-20°C
Size Marker (15-622bp, C109200-100)	-20°C
1x Separation Buffer	*RT
1x Dilution Buffer	*RT
Mineral Oil	*RT

** RT: Room Temperature 30 °C ≥

C. Sample Preparation

- Mix Quantitative Marker (C109109-500Q) and Size Marker (C109200).

Materials	Volume (µl)
Quantitative Marker (C109109-500Q)	5 µl
Size Marker (15-622bp) (C109200)	10 µl
1x Dilution buffer (C104402)	5 µl
Total	20 µl

NOTE: Size Marker mixture should be loaded at regular 0.2ml PCR tube.

- Mix Quantification Marker (C109109-500Q) and Samples.

Materials	Volume (µl)
Quantitative Marker (C109109-500Q)	5 µl
Sample (0.2 ng/ µl to 500 ng/ µl)	X (2~15) µl
1x Dilution buffer (C104402)	15-X (0~13) µl
Total	20 µl



NOTE: Avoid DNA input samples concentration over than specified limits.

- Place size marker mixture at SM/ MA2/ S4 position.



Qsep100


Qsep1




NOTE: Mix all mixtures slightly and spin down before analysis.

D. Software Operation

- Place sample and select corresponding position.
- Enter Sample information (optional).
- Select "Quantification" and enter sample volume (X) at Method Selector.



- Select analytic method "M4-10-06-300" or "M4-10-08-160" in Method Selector.
- Click Run to start analysis.



Contact Information:

Company Name: BiOptic Inc.

Address: (23141) No.108-3, Minquan Rd., Hsin-Tien District, New Taipei City, Taiwan (R.O.C)

Tel: +886-2-2218-8726, Fax: +886-2-2218-8727, E-mail: service@bioptic.com.tw