



Application

Comparative evaluation of reverse transcriptase products for expression analysis of target genes in cultured cells

Product

FastGene® Scriptase II ReadyMix (5X) (LS64)

Manufacturer

NIPPON Genetics EUROPE GmbH

The following data has been provided by the courtesy of Kyoto University customers.

Overview

The comparative study was conducted between FastGene® Scriptase ReadyMix (LS64) and a “two-step RT-qPCR” competitor product, which was currently used. This study was carried out for expression analysis of target genes in cultured cells.

Reverse transcription reaction conditions were evaluated using qPCR as a downstream application, using the current program as it is and obtained cDNA.

Method (reverse transcription reaction)

- RNA sample: RNA recovered from MEF cells (mouse fetal fibroblast), or the same RNA diluted 1:10000 in water
- Reverse transcription kit: ① FastGene® Scriptase II ReadyMix (includes random primer)
② Th cDNA synthesis kit (current product) ※Two-step RT- qPCR kit (including random primer)

• Reaction composition:

① FastGene® Scriptase II ReadyMix

Sample RNA	14 μL
FastGene® Scriptase II ReadyMix	4 μL
dH ₂ O	2 μL ※
total	20 μL

※In this study, dH₂O was added separately, because the same sample RNA solution was used for reaction composition of the current kit ②.

② Th cDNA synthesis kit (current product)

Sample RNA	14 μL
5 x reaction mix	4 μL
10 x enzyme mix	2 μL
total	20 μL

• Reaction conditions

25°C	10 min
42°C	60 min
85°C	5 min
4°C	hold

• Downstream application:

qPCR using SYBR Green I master mix (reagent)
Device used: StepOnePlus™(Thermo Fisher Scientific)

※ For both kits, the reaction conditions of Th's cDNA synthesis kit (current product) were used.

FastGene® Scriptase II Series



The FastGene® Scriptase II reverse transcriptase has a mutation in the reverse transcriptase MMuLV to “improve thermal stability” and “repress the RNase H activity”. This allows more complex applications such as long chain cDNA synthesis for cloning purpose, RT-qPCR, NGS etc.

You can choose from the following three types:

- FastGene® Scriptase II Reverse Transcriptase (LS53)
- FastGene® Scriptase II cDNA Synthesis kit (LS63),
- FastGene® Scriptase II ReadyMix (LS64) Note: For RT-qPCR

FastGene® Scriptase II ReadyMix (5X) (LS64)


This is a “reverse transcription reaction ReadyMix including random primers” released for “quantification of gene expression by qPCR”.

In particular, “quantitative expression of genes by qPCR” often requires processing multiple samples, since this product is a ReadyMix type, it can reduce the time and effort of multiple-sample dispensing operation.





Result



	Sample	Gene expression level	RT kit	
—	1	High	Scriptase II	
	2		company Th kit	
—	1		Middle	Scriptase II
	2			company Th kit
—	1	Low		Scriptase II
	2			company Th kit
—	1		Middle 1:10000 dilution	Scriptase II
	2			company Th kit
—	1	Middle 1:10000 dilution		Scriptase II
	2			company Th kit

Gene expression level	RT kit	Sample	Cr	Cr Mean	Cr SD	Tm1
High	Scriptase II	1	8.35	8.36	0.01	85.26
		2	8.37			85.26
	company Th kit	1	8.54	8.45		85.26
		2	8.37			85.26
Middle	Scriptase II	1	20.14	20.11	0.05	86.45
		2	20.07			86.45
	company Th kit	1	19.81	19.84		86.60
		2	19.86			86.45
Low	Scriptase II	1	28.38	28.40	0.03	83.49
		2	28.42			83.49
	company Th kit	1	29.90	29.75		83.49
		2	29.60			83.49
Middle 1:1000 dilution	Scriptase II	1	34.66	34.22	0.63	86.60
		2	33.78			86.45
	company Th kit	1	33.71	35.36		76.47
		2	37.01			80.79

• Description of target genes:

High : High expression gene

Middle : moderately expression gene

Low : Low expression gene

• Each qPCR measured by duplicates

- Both kits used the reverse transcription reaction program of the current Th company kit. Reverse transcription was performed under the same conditions, and the results of qPCR using the cDNA were compared. As a result, almost identical amplification curves were obtained for both kits, so it was judged that the reverse transcription reaction was stably performed even when the FastGene® Scriptase II ReadyMix was used under the conditions of the current product.
- The reverse transcription of Middle (diluted solution) at a dilution of 1:10000 is expected to be detected about 13.3 cycles later than that of the undiluted solution, but it was clearly confirmed as an experimental result in FastGene® Scriptase II ReadyMix.



Customer's comment

The product was very easy to handle, as I did not use ReadyMix products until now.

In addition, the reverse transcription efficiency of the low expression gene RNA was higher than that of the products conventionally used.