

FastGene® Scriptase Basic

Cat. No.	Pack Size
L\$52	20,000 U
LS52s	2,000 U

Storage:

Store at -20°C.

Description:

FastGene[®] Scriptase Basic, an enhanced MuLV reverse transcriptase, is thermostable at 37-42°C and inactivated at 70°C for 10-15 minutes. It efficiently synthesizes human ACTB, GAPDH, and TFRC cDNAs at 37°C, with optimal temperatures varying by template. It can produce cDNAs up to 5 kb from RNA or DNA templates and has weaker RNase H activity compared to avian myeloblastosis virus RT.



Characteristics:

- Molecular weight: 71 kDa
- Reaction temperature: 37°C or 42°C
- Heat inactivation: 70°C, 10 min

Applications:

- Synthesis of first-strand cDNA
- Array labeling
- cDNA library construction
- 3' and 5' RACE, RT-PCR
- Primer extension

Quality control:

- Purity: >99% on SDS-PAGE
- Endonuclease- and exonuclease-free
- RNase-free
- Inhibitor-free
- Satisfactory cDNA product yield and length



Protocol:

1. Prepare the PCR mixture on ice with the following **components included in this kit**:

Sterile water (RNase free)	up to 20 µ
RNase Inhibitor (40 units/ μ l) - optional	0.5 μl
Primer**	1 μΙ
Template RNA* (not included)	x μl
dNTP mixture (2 mM each)	2 μΙ
FastGene® Scriptase Basic (200 units/µl)	1 µl
10x FastGene [®] Scriptase Basic Buffer	2 µl

* Concentration of template RNA

Total RNA: 1 ng-5 μg Messenger RNA (mRNA): 1 ng-250 ng Specifi RNA: 0.01 ng-0.5 μg ** Prepare one of the following primers Oligo (dT)_{μi}: 50-100 μM Random hexamer: 50-100 μM Specific primer: 15-20 pmol



An additional annealing step is recommended if using random hexamer. Incubate at 25° for 10 min.

2. Incubate at 37°C or 42°C for 60 min. 3. Incubate at 95°C for 5 min to incactivate the reaction.

Note:

FastGene[®] Scriptase Basic can perform optionally over the full range of 42°C-50°C after optimization.

NIPPON Genetics EUROPE GmbH

Mariaweilerstr. 28 a; 52349 Düren; Germany fax: +49 2421 55 496 11, phone: +49 2421 55 496 0 e-mail: info@nippongenetics.eu www.nippongenetics.eu