

FastGene® Scriptase III cDNA 5x ReadyMix

Cat. No.	Pack Size
LS67	100 reactions
LS67s	10 reactions

Storage:

Store at -20°C.

Description:

FastGene® Scriptase III cDNA 5x ReadyMix is an improved version of M-MLV Reverse Transcriptase, which exhibits high thermal stability (~55°C) and low RNase H activity. It is a mixture of FastGene® Scriptase III, RNase inhibitor, MgCl₂, and dNTPs. It is ready to use by adding template RNA and primers. The synthesis of cDNA can be carried out using either total RNA or poly(A)⁺-selected RNA, with individual priming options including oligo(dT), random primers, or a gene-specific primer.

Characteristics:

- Produces excellent yields of amplifiable cDNA.
- Can synthesize long cDNA targets, up to 12 kb, from small amounts of input RNA.
- Able to reverse transcribe through RNA secondary structures.
- Strict enzyme purity specifications.

Applications:

- First-Strand cDNA synthesis
- Conventional PCR
- Real-time quantitative RT-PCR (qRT-PCR)

Quality control:

- Purity: >99% on SDS-PAGE
- Endonuclease- and exonuclease-free
- RNase-free
- Inhibitor-free
- Satisfactory yield and length of cDNA products

Protocol:

1. Prepare the PCR mixture on ice with the following **components included in this kit**:

FastGene® Scriptase III cDNA 5x ReadyMix	4 µl
Primer (80 µM oligo(dT), or 2 µM gene-specific primer, or 50 ng/µL random hexamers)	1 µl
Template RNA* (not included)	x µl
Sterile water (RNase free)	up to 20 µl

2. Incubate at 25°C for 10 min and 42°C for 60 min.

3. Incubate at 85° for 5 min to inactivate the reaction.

Note:

- High quality RNA is needed for accurate quantification in qPCR. RNA should be kept in the absence of RNase contamination.
- Template RNA can range up to 2.5 µg in a 20 µL cDNA synthesis reaction.
- Amplification grade DNase I can be used to remove genomic DNA contamination from the total RNA.
- A shorter incubation time of minimum 30 min and higher synthesis temperatures (up to 50°C) can be used without a loss of performance.
- Longer incubation times may be used for increased yields of cDNA.
- Undiluted synthesized cDNA may be used up to 10% of the qPCR reaction volume (e.g., for a 20 µL qPCR, use up to 2 µL of undiluted cDNA).

NIPPON Genetics EUROPE GmbH

Mariaweyerstr. 28 a; 52349 Düren; Germany
fax: +49 2421 55 496 11, phone: +49 2421 55 496 0
e-mail: info@nippongenetics.eu
www.nippongenetics.eu