

## Compatibility list of the FastGene® PCR Plates

### FastGene® 96-well Plate

- Non-skirted
- Volume: 0.2 ml
- **Cat. No.: FG-170225 (25 Plates)**



Supplier	Type
NIPPON Genetics EUROPE	FastGene® Ultra Cycler Gradient (FG-TC01) FastGene® Ultra Cycler Uniform (FG-TC02)
Applied Biosystems (ABI)	GeneAmp™ 2720 GeneAmp™ 9600 GeneAmp™ 9700 7000, 7300, 7500, 7700
Biometra	T1 Thermocycler Tgradient
Bio-Rad	PTC-200 DNA Engine™ PTC-225 DNA Tetrad
Gene Technologies	G-Storm GS1
Stratagene	M x 4000™ M x 3000™
Takara	Dice

### FastGene® 96-well Plate

- Semi-skirted
- Volume: 0.3 ml
- **Cat. No.: FG-190250 (50 Plates)**



Supplier	Type
Applied Biosystems (ABI)	GeneAmp™ 2700/2720 GeneAmp™ 9600 GenAmp™ 9700 7000/7300/7500/7700/7900
Biometra	T1 Thermocycler Tgradient Tadvanced
Bio-Rad	iCycler™ iQ / iQ4 / iQ5 MyCycler C1000 S1000 PTC100™ PTC200 (DNA Engine™) PTC220/221 (DNA Dyad™) PTC225 (DNA Tetrad™) MyiQ
Eppendorf	Mastercycler™ Mastercycler™ EP Gradient Mastercycler™ ep realplex
MWG	The Q-Lifecycler
PeqLab	peqSTAR 96
Technne	Touchgene Gradient Flexigene Genius

## FastGene® Fast 96-well Plate

- Semi-skirted
- Volume: 0.1 ml
- **Cat. No.: FG-03890 (25 Plates)**



Supplier	Type
NIPPON Genetics EUROPE	FastGene® Ultra Cycler Gradient (FG-TC01) FastGene® Ultra Cycler Uniform (FG-TC02)
Applied Biosystems (ABI)	StepOnePlus™ ABI 7500 Fast ABI 9800 QuantStudio™

## FastGene® Fast 96-well Plate

- Full-skirted
- Volume: 0.2 ml
- **Cat. No.: FG-180230 (30 Plates)**



Supplier	Type
NIPPON Genetics EUROPE	FastGene® Ultra Cycler Gradient (FG-TC01) FastGene® Ultra Cycler Uniform (FG-TC02)
Biometra	Uno T1 Thermocycler Tgradient Trobot Tadvanced
Bio-Rad	C1000 S1000 PTC100™ PTC200 (DNA Engine™) PTC220/221 (DNA Dyad™) PTC225 (DNA Tetrad™) CFX96 Opticon™ Opticon2™ Chromo4™
Corbett Research	Palm Cyclor
Eppendorf	Master-Cycler® Gradient Master-Cycler® EP Gradient Mastercycler™ ep realplex
G-Storm	GS1, GS2, GS3
MWG	Primus 96 The Q-Lifecycler
Stratagene	Gradient Cycler
Techne	TC412/512 Touchgene Gradient Flexigene Genius
Thermo Fisher Scientific	PCR Express/Px2/PxE Multiblock System MBS Touchdown Omnigene Omni-E

## FastGene® Fast 96-well Plate Roche® Style

- Semi-skirted
- Volume: 0.2 ml
- **Cat. No.: FG-210250 (50 Plates)**



Supplier	Type
Analytik-Jena	qTOWER 2.0
Applied Biosystems (ABI)	StepOne Plus
Bio-Rad	Opticon Opticon2 Chromo4 CFX96 CFX 96 touch
Roche	LightCycler 480 LightCycler 96

## FastGene® Fast 96-well Plate FROSTED ABI® Style

- Semi-skirted
- Volume: 0.2 ml
- **Cat. No.: FG-200250 (50 Plates)**



Supplier	Type
NIPPON Genetics EUROPE	FastGene® Ultra Cyclor Gradient (FG-TC01) FastGene® Ultra Cyclor Uniform (FG-TC02)
Applied Biosystems (ABI)	7000, 7300, 7500, 7700, 7900 7900 HAT QuantStudio™ ViiA7™ 310, 3100 Genetic Analyser 3130 Genetic Analyser 3500XL 3700, 3730, 3730XL DNA Analyser GeneAmp® 2700, 2720, 9600 GeneAmp® 9700 Veriti 0.2ml (96 well block)
Eppendorf	Mastercycler™ ep realplex
Statagene/Agilent	Mx4000™

## FastGene® Fast 384-well Plate

- Full-skirted
- Volume: 0.05 ml
- **Cat. No.: FG-300150 (50 Plates)**



Supplier	Type
Agilent	SureCycler 8800
Amersham	MegaBACE 4000
Analytik Jena	FlexCycler T1 Thermal Cycler Tprofessional/Basic Tprofessional/Standard Tprofessional/TRIO Trobot Uno II
Applied Biosystems (ABI)	9700, Veriti/Fast 3100 3130 3700 3730/3730x 7900HT/Fast StepOne Plus Via7/Fast
Bio-Rad	C1000 DNA Engine Dyad Dyad Disciple Engine Tetrad 2 PTC-100 PTC-200 PTC-225 Tetrad S1000 T100 CFX384 Touch
Corbett	Palm-Cycler 384
Eppendorf	Mastercycler M384
MWG	Primus 384
PeqLab	PeqStar 384
Qiagen	FlexiGene
SensoQuest	Thermocycler 384
Techne	Genius Genius Quad Prime Prime Elite Prime Elite Satellite PrimeG Touchgene Gradient Touchgene X
Thermo Hybrid	Multiblock System Omnigene PCR Express Px2 PxE Touchdown